

# SYNERGY™

CT PICC



Simple grip luers are easy on gloves and maintain strength and crack resistance

Polyurethane catheter construction provides radiopacity, durability and strength

Easy legible identification on clamps and extensions, indicating power injectability and maximum pressure

Low profile hub allows for suture or sutureless securement



SYNERGY™  
CT PICC

THE SYNERGY™ CT PICC OFFERS THE EASE AND EFFECTIVENESS OF PICC ACCESS WITH POWER INJECTABILITY, CREATING A SYNERGY OF INTRAVENOUS THERAPY OPTIONS FOR HEALTHCARE PRACTITIONERS.

# SYNERGY™ CT PICC

## SYNERGY™ CT Peripherally Inserted Central Catheters (PICC)

Catalogue Number	PICC Description	Guide Wire	Introducer	Box Quantity
<b>SYNERGY™ CT Peripherally Inserted Central Catheters: Nurses Kit</b>				
A14-04160-N	4f x 60cm Single Lumen SYNERGY™ CT PICC	45cm Nitinol GW W/Tungsten Tip	4.5f x 05cm PTFE Tearaway	5
A14-05160-N	5f x 60cm Single Lumen SYNERGY™ CT PICC	45cm Nitinol GW W/Tungsten Tip	5.5f x 05cm PTFE Tearaway	5
A14-05255-N	5f x 55cm Dual Lumen SYNERGY™ CT PICC	45cm Nitinol GW W/Tungsten Tip	5.5f x 05cm PTFE Tearaway	5
A14-05260-N	5f x 60cm Dual Lumen SYNERGY™ CT PICC	45cm Nitinol GW W/Tungsten Tip	5.5f x 05cm PTFE Tearaway	5
A14-06260-N	6f x 60cm Dual Lumen SYNERGY™ CT PICC	45cm Nitinol GW W/Tungsten Tip	6.5f x 05cm PTFE Tearaway	5
A14-06360-N	6f x 60cm Triple Lumen SYNERGY™ CT PICC	45cm Nitinol GW W/Tungsten Tip	6.5f x 05cm PTFE Tearaway	5
<b>SYNERGY™ Peripherally Inserted Central Catheters: Basic Kit</b>				
A14-04150	4f x 50cm Single Lumen SYNERGY™ CT PICC	70cm Nitinol GW W/Tungsten Tip	4.5f x 10cm PTFE Tearaway	5
A14-05155	5f x 55cm Single Lumen SYNERGY™ CT PICC	70cm Nitinol GW W/Tungsten Tip	5.5f x 10cm PTFE Tearaway	5
A14-05255	5f x 55cm Dual Lumen SYNERGY™ CT PICC	70cm Nitinol GW W/Tungsten Tip	5.5f x 10cm PTFE Tearaway	5
A14-06260	6f x 60cm Dual Lumen SYNERGY™ CT PICC	70cm Nitinol GW W/Tungsten Tip	6.5f x 10cm PTFE Tearaway	5
A14-06360	6f x 60cm Triple Lumen SYNERGY™ CT PICC	70cm Nitinol GW W/Tungsten Tip	6.5f x 10cm PTFE Tearaway	5
<b>SYNERGY™ Peripherally Inserted Central Catheters: IR Kit</b>				
A14-04160-IR	4f x 60cm Single Lumen SYNERGY™ CT PICC	130cm Nitinol GW W/Tungsten Tip	4.5f x 10cm PTFE Tearaway	5
A14-05260-IR	5f x 60cm Dual Lumen SYNERGY™ CT PICC	130cm Nitinol GW W/Tungsten Tip	5.5f x 10cm PTFE Tearaway	5
A14-06260-IR	6f x 60cm Dual Lumen SYNERGY™ CT PICC	130cm Nitinol GW W/Tungsten Tip	6.5f x 10cm PTFE Tearaway	5
A14-06360-IR	6f x 60cm Triple Lumen SYNERGY™ CT PICC	130cm Nitinol GW W/Tungsten Tip	6.5f x 10cm PTFE Tearaway	5

### All Kits Contain The Following:

SYNERGY™ CT PICC	Safety Scapel
Needle Free Valve's	24" Paper Tape Measure
0.15 x 75cm Twisted Wire Stylet	PICC Stylet TAG
21GA x 7cm Safety Introducer Needle w/ Eco Tip	HLIC PICC Patient Info Guide
PTFE Tearaway Introducer	HLIC CT PICC Chart Sticker
Nitinol Straight Mandrel GW w/ Tungsten Tip	HLIC CT PICC IFU
10cc Luer Lock Syringe	Statlock® Securement Device w/ Fixed or Sliding POS

Premium Nitinol Guide Wires and Eco Tip Safety Introducer Needles in all kits

The SYNERGY™ CT PICC is available in 4F single lumen through 6F triple lumen catheter configurations. Nursing, Basic, and IR Insertion trays available for all PICC sizes.



**Manufactured By:**  
Health Line International Corp.  
803 N. 1250 W. #1  
Centerville, Utah  
USA / 84014  
Phone: 1.877.847.4542  
www.hlic.net

**European Authorized Rep.:**  
Obelis s.a.  
Boulevard General Wahis 53  
1030 Brussels, Belgium  
Phone: +(32)2.732.59.54  
Fax: +(32)2.732.60.03

**CE: 0197**