


AQUEOUS™

Peritoneal Dialysis Catheters and Kits



Health Line's AQUEOUS™ Peritoneal Dialysis (PD) catheters are used to treat renal failure and help patients get back on track to living a healthy life. They are indicated for use in long term (more than 30 days) therapy. They're available in a variety of sizes to accommodate neonatal, pediatric, and adult patients. The catheters are offered in coiled, straight, and swan neck designs so clinics have options necessary to meet all the patients needs.



STRATEGICALLY PLACED SIDE HOLES
allow for increased flow and reduced
omental occlusions

POLYESTER CUFF
allows for secure tunneled anchoring

Features and Benefits

Radiopaque Silicone - Provides excellent visibility and biocompatibility.

Standard Trays and Components - Trays are organized with familiar industry components for clinical safety and user friendly procedural steps; including spreadable sheath introducer, double ended guidewire, tunneling stylet, and introducer needle.


Wide Range of Use - PD Catheters are designed to allow easy use in the home, at work, or even during travel, and various size and shape options ensure usability in Neonatal, Pediatric, and Adult patients in any situation or need.

Aqueous™ Peritoneal Dialysis Catheters and Kits

Description	Catalogue No.	Box Qty.	Carton Qty.	French/mm	Length (cm)	No. of Cuffs	Description
Peritoneal Straight Catheters (Adult, Pediatric and Neonatal)	PCS-42-1	5	20	14.7/4.9	42	1	Straight
	PCS-42-2	5	20	14.7/4.9	42	2	Straight
	PCS-47-2	5	20	14.7/4.9	47	2	Straight
	PCS-31-1N	5	20	14.7/4.9	31	1	Straight Neonatal
	PCS-32-2P	5	20	14.7/4.9	32	2	Straight Pediatric
	PCS-37-2P	5	20	14.7/4.9	37	2	Straight Pediatric
Peritoneal Coiled Catheters (Adult and Pediatric)	PCC-39.25-1P	5	20	14.7/4.9	39	1	Coiled Pediatric
	PCC-56.5-2	5	20	14.7/4.9	57	2	Coiled
	PCC-62-2	5	20	14.7/4.9	62	2	Coiled
Peritoneal Swan Neck Catheters	PCC-43-2	5	20	14.55/4.85	43	2	Swan Neck
	PCC-37.5-2P	5	20	14.55/4.85	37.5	2	Swan Neck Pediatric
Peritoneal Placement Kits	PC-K	5	20				

Package Contents

Component	Quantity
Peritoneal Silicone Catheter	1
Female luer	1
Male luer cap	1
Catheter clip	1



PC-K Placement Kit

Peritoneal Catheter Placement Kit Contents

Component	Quantity	Component	Quantity
Introducer needle, thin-wall, 18 gauge x 2 3/4" (7cm) length	1	Spreadable sheath introducer, 16fr	1
		Tunneling stylet	1
"J" Double-ended guidewire, 0.038" O.D. x 27.5" (70cm) length	1	#11 Scalpel	1
		Gauze sponges	6
Dilator, 10 French (3.3mm) O.D. x 6" (15.2cm) length	1	Syringe, 10cc	1
		Female luer	1
Catheter clip	1	Male luer cap	1

For ordering information or for questions about the Aqueous™, please contact your Health Line Sales Consultant at (877) 847-4542 and visit us at www.hlic.net.





Health Line International Corp.
5675 West 300 South
Salt Lake City, Utah 84104

Phone: (877) 847-4542
Fax: (855) 228-1336
www.hlic.net

ISO13485
Certified Quality System
MKT-100-07 Rev. 01